


PRODUCT DATA SHEET

Revision 1	384 Well Cell Culture Plate, Black with Clear F- bottom	
Validfor Item-No.:	1015700/1015604/1015606	

1.	Description / Specification	
1.1	Description	Cell Culture Plate, 384 well, black with clear flat bottom , with clear lid , alphanumeric well coding, single position lid with condensation rings, sterile 1015700 384 well, Cell culture plate, Non-TC treated, sterile 1015604 384 well, Cell culture plate, TC treated, sterile 1015606 384 well, Cell culture plate, TC treated, sterile
1.2	Dimensions	Plate : see Customer Drawing 1015700/1015604:Individually Wrapped in Peelable 1015606:Individually paper-plastic wrapped
1.3	Volume	Max volume: 0.112 ml (mathematical calculated) Working volume:0.02-0.08 ml Growth area / well: 0.06 cm ²
1.4	Material / Resin	Plate : PS (Polystyrene) Lid: PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	Plate : black with clear bottom Lid: clear
1.6	Sterilization	SAL 10 ⁻⁶
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	50
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
4.1	Laboratory use only	

Prior Issue	Drawn	Approved	Released	
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.