


PRODUCT DATA SHEET

Revision 1	Laminin-Coated Surface, Cell Culture Plate	
Validfor Item-No.:	1130002/1132002/1134002	

1.	Description / Specification	
1.1	Description	Laminin-Coated surface, Cell Culture Plate, Sterile 1130002 Laminin-Coated surface, cell culture plate, 6 well, transparent flat bottom, sterile 1132002 Laminin-Coated surface, cell culture plate, 24 well, transparent flat bottom, sterile 1134002 Laminin-Coated surface, cell culture plate, 96 well, transparent flat bottom, sterile
1.2	Dimensions	<u>Plate</u> : see Customer Drawing
1.3	Volume	Max. volume: 4ml (6 well) / 1ml (24 well) / 0.2ml (96 well) Working volume: 4ml (6 well) / 1ml (24 well) / 0.1ml (96 well) Growth area: 9.6 cm ² (6 well) / 1.9 cm ² (24 well) / 0.32 cm ² (96 well)
1.4	Material / Resin	<u>Plate</u> : PS (Polystyrene) <u>Lid</u> : PS (Polystyrene) The materials for manufacturing are <i>USP Class VI</i>
1.5	Colour	<u>Plate</u> : clear <u>Lid</u> : clear
1.6	Sterilization	SAL 10 ⁻⁶
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. Printed on Packaging

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: +2 °C to +8 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	5
3.2	Pieces / Box	5
3.3	Lot-No.	DDD YY XXX (day, year, product No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Laboratory use only

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.