


PRODUCT DATA SHEET

Revision 1/	50µl Automation Tips (Adaptation of Tecan)	
Validfor Item-No.:	3020001/3020011/3020101/3020111/3020002/ 3020012/3020102/3020112	

1.	Description / Specification	
1.1	Description	Automation Tips (Adaptation of Tecan), 50µl, clear, sterile, boxed 3020001 Automation Tips (Adaptation of Tecan), 50µl, universal, clear, sterile, boxed 3020011 Automation Tips (Adaptation of Tecan), 50µl, universal, clear, filter, sterile, boxed 3020101 Automation Tips (Adaptation of Tecan), 50µl, low retention, clear, sterile, boxed 3020111 Automation Tips (Adaptation of Tecan), 50µl, low retention, clear, filter, sterile, boxed 3020002 Automation Tips (Adaptation of Tecan), 50µl, universal, clear, sterile, boxed 3020012 Automation Tips (Adaptation of Tecan), 50µl, universal, clear, filter, sterile, boxed 3020102 Automation Tips (Adaptation of Tecan), 50µl, low retention, clear, sterile, boxed 3020112 Automation Tips (Adaptation of Tecan), 50µl, low retention, clear, filter, sterile, boxed
1.2	Dimensions	<u>Flask</u> : see Customer Drawing
1.3	Volume	See 1.1Description
1.4	Material / Resin	Pipette tip: PP (Polypropylene), free of heavy metal Filter: PE (Polyethylene) Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: clear Filter: white Rack/lid: clear
1.6	Sterilization	SAL 10 ⁻⁶
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only Expiry date and Lot-No. printed on the outer box

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +121 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	3 years
2.6	Other Information	-

3.	Packaging	3020001/3020011/3020101/3020111	3020002/3020012/3020102/3020112
3.1	Pieces / Box	96	96
3.2	Pieces / middle Box	10	5
3.3	Pieces / Case	5	20
3.4	Lot-No.	DDD YY XXX (day, year, product No.)	
3.5	Other Information	Certificate of Quality to download	

4.	Other Information	
4.1	Laboratory use only	

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Zhejiang Saining Biotechnology Co., Ltd. This document may not be reproduced for any reason without written permission from Zhejiang Saining Biotechnology Co., Ltd. All rights of design, invention, and copyright are reserved.
Revision 1	Date 1 June 2024	Date 1 June 2024	Date 1 June 2024	
Date 01.06.2024	Name Chen Anjing	Name Zhang Heng	Name Wang Kedong	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.